

[KO] CD147 Mouse mAb

ABZ1187

Product Overview

Name	[KO] CD147 Mouse mAb
Accession(Primary)	P35613
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human BSG produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 55 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200

Target information (P35613)

Gene Name	BSG
About Gene Symbol	BSG (Basigin (Ok Blood Group)) is a collagen extracellular matrix protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.