

[KO] TAK1 Mouse mAb

ABZ1188

Product Overview

Name	[KO] TAK1 Mouse mAb
Accession(Primary)	O43318
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human MAP3K7 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 67 kDa; Observed MW: 80 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (O43318)

Gene Name	MAP3K7
About Gene Symbol	MAP3K7 (Mitogen-Activated Protein Kinase Kinase 7) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and TGF- β signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.