

[KO] TARBP2 Mouse mAb

ABZ1189

Product Overview

Name	[KO] TARBP2 Mouse mAb
Accession(Primary)	Q15633
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human TARBP2 produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 50 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q15633)

Gene Name	TARBP2
About Gene Symbol	TARBP2 (TARBP2 Subunit Of RISC Loading Complex) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and i