

[KO] TDG Mouse mAb

ABZ1190

Product Overview

Name	[KO] TDG Mouse mAb
Accession(Primary)	Q13569
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human TDG protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 55 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q13569)

Gene Name	TDG
About Gene Symbol	