

[KO] Thimet Oligopeptidase Mouse mAb

ABZ1191

Product Overview

Name	[KO] Thimet Oligopeptidase Mouse mAb
Accession(Primary)	P52888
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human THOP1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 79 kDa; Observed MW: 79 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P52888)

Gene Name	THOP1
About Gene Symbol	