

[KO] TSN Mouse mAb

ABZ1192

Product Overview

Name	[KO] TSN Mouse mAb
Accession(Primary)	Q15631
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human TSN protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q15631)

Gene Name	TSN
About Gene Symbol	