

[KO] Vimentin Mouse mAb

ABZ1193

Product Overview

Name	[KO] Vimentin Mouse mAb
Accession(Primary)	P08670
Host Species	Mouse
Isotype	IgG2a
Immunogen	Purified recombinant human VIM protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P08670)

Gene Name	VIM
About Gene Symbol	VIM (Vimentin) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.