

[KO] USB1 Mouse mAb

ABZ1194

Product Overview

Name	[KO] USB1 Mouse mAb
Accession(Primary)	Q9BQ65
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human USB1 produced in E.coli.
Calculated MW	Calculated MW: 30 kDa; Observed MW: 34 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q9BQ65)

Gene Name	USB1
About Gene Symbol	