

[KO] XLF Mouse mAb

ABZ1195

Product Overview

Name	[KO] XLF Mouse mAb
Accession(Primary)	Q9H9Q4
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human NHEJ1 produced in 293T cell.
Calculated MW	Calculated MW: 33 kDa; Observed MW: 37 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (Q9H9Q4)

Gene Name	NHEJ1
About Gene Symbol	NHEJ1 (Non-Homologous End Joining Factor 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.