

[KO] TXNDC5 Mouse mAb

ABZ1198

Product Overview

Name	[KO] TXNDC5 Mouse mAb
Accession(Primary)	Q8NBS9
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human TXNDC5 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 55 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q8NBS9)

Gene Name	TXNDC5
About Gene Symbol	TXNDC5 (Thioredoxin Domain Containing 5) is a protein kinase and signaling regulator involved in angiogenesis. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechan