

[KO] CLPP Mouse mAb

ABZ1199

Product Overview

Name	[KO] CLPP Mouse mAb
Accession(Primary)	Q16740
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human CLPP produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 25 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q16740)

Gene Name	CLPP
About Gene Symbol	CLPP (Caseinolytic Mitochondrial Matrix Peptidase Proteolytic Subunit) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate