

[KO] USP9X Mouse mAb

ABZ1200

Product Overview

Name	[KO] USP9X Mouse mAb
Accession(Primary)	Q93008
Host Species	Mouse
Isotype	IgG2a
Immunogen	Purified recombinant human USP9X protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 290 kDa; Observed MW: 290 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q93008)

Gene Name	USP9X
About Gene Symbol	USP9X (Ubiquitin Specific Peptidase 9 X-Linked) is a deubiquitinase involved in ubiquitin-proteasome pathway. It is an important research target in neurodegenerative disease and is widely studied to elucidate cellular mechanisms and identify therapeutic t