

[KO] ASCC2 Mouse mAb

ABZ1202

Product Overview

Name	[KO] ASCC2 Mouse mAb
Accession(Primary)	Q9H1I8
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human ASCC2 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 86 kDa; Observed MW: 110 kDa
Recommended Dilution	WB: 1/500-1/1000; FC: 1/50-1/100

Target information (Q9H1I8)

Gene Name	ASCC2
About Gene Symbol	ASCC2 (Activating Signal CointegRator 1 Complex Subunit 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.