

[KO] PPM1D Mouse mAb

ABZ1203

Product Overview

Name	[KO] PPM1D Mouse mAb
Accession(Primary)	O15297
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human PPM1D protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 67 kDa; Observed MW: 80 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (O15297)

Gene Name	PPM1D
About Gene Symbol	