

[KO] PRAK Mouse mAb

ABZ1204

Product Overview

Name	[KO] PRAK Mouse mAb
Accession(Primary)	Q8IW41
Host Species	Mouse
Isotype	IgG2a
Immunogen	Purified recombinant human MAPKAPK5 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 54 kDa; Observed MW: 54 kDa
Recommended Dilution	WB: 1/500-1/1000

Target information (Q8IW41)

Gene Name	MAPKAPK5
About Gene Symbol	