

[KO] PDXK Mouse mAb

ABZ1205

Product Overview

Name	[KO] PDXK Mouse mAb
Accession(Primary)	O00764
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human PDXK produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (O00764)

Gene Name	PDXK
About Gene Symbol	