

[KO] Beclin 1 Mouse mAb

ABZ1206

Product Overview

Name	[KO] Beclin 1 Mouse mAb
Accession(Primary)	Q14457
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human BECN1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 65 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (Q14457)

Gene Name	BECN1
About Gene Symbol	BECN1 (Beclin 1) is an autophagy regulatory protein involved in autophagy and apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.