

[KO] SRSF9 Mouse mAb

ABZ1207

Product Overview

Name	[KO] SRSF9 Mouse mAb
Accession(Primary)	Q13242
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human SRSF9 produced in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q13242)

Gene Name	SRSF9
About Gene Symbol	SRSF9 (Serine And Arginine Rich Splicing Factor 9) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.