

[KO] alpha Actinin Mouse mAb

ABZ1208

Product Overview

Name	[KO] alpha Actinin Mouse mAb
Accession(Primary)	P12814
Host Species	Mouse
Isotype	IgG2a
Immunogen	Purified recombinant human ACTN1 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 103 kDa; Observed MW: 103 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; ICC: 1/100-1/200

Target information (P12814)

Gene Name	ACTN1
About Gene Symbol	