

[KO] CBR1 Mouse mAb

ABZ1209

Product Overview

Name	[KO] CBR1 Mouse mAb
Accession(Primary)	P16152
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human CBR1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 30 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P16152)

Gene Name	CBR1
About Gene Symbol	