

## [KO] SELENBP1 Mouse mAb

ABZ1211

### Product Overview

|                             |                                                                                |
|-----------------------------|--------------------------------------------------------------------------------|
| <b>Name</b>                 | [KO] SELENBP1 Mouse mAb                                                        |
| <b>Accession(Primary)</b>   | Q13228                                                                         |
| <b>Host Species</b>         | Mouse                                                                          |
| <b>Isotype</b>              | IgG2a                                                                          |
| <b>Immunogen</b>            | Full length human recombinant protein of human SELENBP1 produced in 293T cell. |
| <b>Species Reactivity</b>   | Rat                                                                            |
| <b>Calculated MW</b>        | Calculated MW: 52 kDa; Observed MW: 52 kDa                                     |
| <b>Recommended Dilution</b> | WB: 1/500-1/1000; IHC: 1/50-1/100                                              |

### Target information (Q13228)

|                          |                                                                                                                                                                                                                                                                                      |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | SELENBP1                                                                                                                                                                                                                                                                             |
| <b>About Gene Symbol</b> | SELENBP1 (Selenium Binding Protein 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |