

[KO] PAPSS2 Mouse mAb

ABZ1212

Product Overview

Name	[KO] PAPSS2 Mouse mAb
Accession(Primary)	O95340
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human PAPSS2 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (O95340)

Gene Name	PAPSS2
About Gene Symbol	