

[KO] USP10 Mouse mAb

ABZ1213

Product Overview

Name	[KO] USP10 Mouse mAb
Accession(Primary)	Q14694
Host Species	Mouse
Isotype	IgG2a
Immunogen	Full length human recombinant protein of human USP10 produced in 293T cell.
Calculated MW	Calculated MW: 87 kDa; Observed MW: 120 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (Q14694)

Gene Name	USP10
About Gene Symbol	