

[KO] TRIP13 Mouse mAb

ABZ1215

Product Overview

Name	[KO] TRIP13 Mouse mAb
Accession(Primary)	Q15645
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant human TRIP13 protein fragments expressed in E.coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 40 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (Q15645)

Gene Name	TRIP13
About Gene Symbol	