

[KO] ERp72 Mouse mAb

ABZ1217

Product Overview

Name	[KO] ERp72 Mouse mAb
Accession(Primary)	P13667
Host Species	Mouse
Isotype	IgG2b
Immunogen	Full length human recombinant protein of human PDIA4 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 73 kDa; Observed MW: 68 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100

Target information (P13667)

Gene Name	PDIA4
About Gene Symbol	PDIA4 (Protein Disulfide Isomerase Family A Member 4) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.