

## [KO] CD68 Mouse mAb

ABZ1218

### Product Overview

|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | [KO] CD68 Mouse mAb                                                       |
| <b>Accession(Primary)</b>   | P34810                                                                    |
| <b>Host Species</b>         | Mouse                                                                     |
| <b>Isotype</b>              | IgG2b                                                                     |
| <b>Immunogen</b>            | Purified recombinant human CD68 protein<br>fragments expressed in E.coli. |
| <b>Calculated MW</b>        | Calculated MW: 37 kDa; Observed MW: 68 kDa                                |
| <b>Recommended Dilution</b> | WB: 1/500-1/1000; IHC: 1/50-1/100                                         |

### Target information (P34810)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | CD68 |
| <b>About Gene Symbol</b> |      |