

[KO] Keap1 Mouse mAb

ABZ1219

Product Overview

Name	[KO] Keap1 Mouse mAb
Accession(Primary)	Q14145
Host Species	Mouse
Isotype	IgG1
Immunogen	Full length human recombinant protein of human KEAP1 produced in 293T cell.
Species Reactivity	Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 70 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; FC: 1/50-1/100

Target information (Q14145)

Gene Name	KEAP1
About Gene Symbol	KEAP1 (Kelch Like ECH Associated Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.