

[KO] MAGEA4 Mouse mAb

ABZ1220

Product Overview

Name	[KO] MAGEA4 Mouse mAb
Accession(Primary)	P43358
Host Species	Mouse
Isotype	IgG1
Immunogen	Purified recombinant fragment of human MAGEA4 expressed in E. Coli.
Species Reactivity	Rat
Calculated MW	Calculated MW: 35kDa;Observed MW: 43kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/100-1/200; FC: 1/50-1/100; ELISA: 1/10000

Target information (P43358)

Gene Name	MAGEA4
About Gene Symbol	