

[KO] MMP9 Rabbit pAb

ABZ1221

Product Overview

Name	[KO] MMP9 Rabbit pAb
Accession(Primary)	P14780
Host Species	Rabbit
Isotype	IgG
Immunogen	A synthetic peptide of human MMP9
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 78 kDa; Observed MW: 100 kDa
Recommended Dilution	WB: 1/500-1/1000; IHC: 1/50-1/100; IF: 1/50-1/200

Target information (P14780)

Gene Name	MMP9
About Gene Symbol	