

MD2 Rabbit pAb

ABZ1273

Product Overview

Name	MD2 Rabbit pAb
Accession(Primary)	Q9Y6Y9
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:5000-1:20000

Target information (Q9Y6Y9)

Gene Name LY96**About Gene Symbol**

