

RANK Rabbit pAb


ABZ1274

Product Overview


Name	RANK Rabbit pAb
Accession(Primary)	Q9Y6Q6
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 66 kDa; Observed MW: 66 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:500-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200

Target information (Q9Y6Q6)

Gene Name	TNFRSF11A
About Gene Symbol	



Western blot analysis of RANK in Hela lysates using RANK antibody.



Immunohistochemistry analysis of paraffin-embedded Human liver using RANK antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.