


B3GALT1 Rabbit pAb

ABZ1275

Product Overview

Name	B3GALT1 Rabbit pAb
Accession(Primary)	Q9Y5Z6
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (Q9Y5Z6)

Gene Name B3GALT1**About Gene Symbol**

Western blot analysis of beta 3 Gal T1 in various lysates using beta 3 Gal T1 antibody.

