

Toll-Like Receptor 6 Rabbit pAb

ABZ1277


Product Overview

Name	Toll-Like Receptor 6 Rabbit pAb
Accession(Primary)	Q9Y2C9
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:200-1:500; IF-1:50-1:200

Target information (Q9Y2C9)

Gene Name TLR6

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded mouse Spleen Tissue using Toll-Like Receptor 6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.