

## CaMKII alpha/delta Rabbit pAb


ABZ1278

### Product Overview


|                             |                                                                     |
|-----------------------------|---------------------------------------------------------------------|
| <b>Name</b>                 | CaMKII alpha/delta Rabbit pAb                                       |
| <b>Accession(Primary)</b>   | Q9UQM7                                                              |
| <b>Host Species</b>         | Rabbit                                                              |
| <b>Clonality</b>            | Polyclonal                                                          |
| <b>Isotype</b>              | IgG                                                                 |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                   |
| <b>Calculated MW</b>        | Calculated MW: 54 kDa; Observed MW: 54 kDa                          |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:200-1:1000; ELISA-1:10000 |

### Target information (Q9UQM7)


|                          |               |
|--------------------------|---------------|
| <b>Gene Name</b>         | CAMK2A/CAMK2D |
| <b>About Gene Symbol</b> |               |




Western blot analysis of CaMKII alpha/delta in NIH3T3/K562 lysates using CaMKII alpha/delta antibody.



Immunofluorescence analysis of CaMKII alpha/delta in MCF-7 using CaMK2 alpha/delta antibody. The picture on the right is blocked using the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded Human brain tissue using CaMK2 alpha/delta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of CaMKII alpha/delta in NIH/3T3 and K562 lysates using CaMKII alpha/delta antibody. The lane on the right is blocked with the synthesized peptide.