

HDAC5 Rabbit pAb


ABZ1279

Product Overview


Name	HDAC5 Rabbit pAb
Accession(Primary)	Q9UQL6
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 122 kDa; Observed MW: 122 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (Q9UQL6)


Gene Name	HDAC5
About Gene Symbol	




Western blot analysis of HDAC5 in various lysates using HDAC5 antibody.



Western blot analysis of HDAC5 in K562 lysates using HDAC5 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue using HDAC5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of HDAC5 in colo lysates using HDAC5 antibody. The lane on the right is blocked with the synthesized peptide.