

ASC Rabbit pAb


ABZ1282

Product Overview

Name	ASC Rabbit pAb
Accession(Primary)	Q9ULZ3
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 22 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:40000

Target information (Q9ULZ3)

Gene Name	PYCARD
About Gene Symbol	



Western blot analysis of ASC in various lysates
using ASC antibody.