

Rab26 Rabbit pAb


ABZ1283

Product Overview


Name	Rab26 Rabbit pAb
Accession(Primary)	Q9ULW5
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:40000

Target information (Q9ULW5)


Gene Name	RAB26
About Gene Symbol	



Western blot analysis of Rab26 in K562 lysates using Rab 26 antibody.



Immunohistochemistry analysis of paraffin-embedded rat brain using Rab26 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded rat brain using Rab26 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.