


C1RL Rabbit pAb

ABZ1285

Product Overview

Name	C1RL Rabbit pAb
Accession(Primary)	Q9NZP8
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (Q9NZP8)

Gene Name C1RL**About Gene Symbol**

Western blot analysis of C1RL in hela lysates using C1RL antibody.

