

TREM2 Rabbit pAb


ABZ1287

Product Overview


Name	TREM2 Rabbit pAb
Accession(Primary)	Q9NZC2
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25-35-50 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:5000-1:20000; IHC-1:50-1:100; IF-1:50-1:200

Target information (Q9NZC2)


Gene Name	TREM2
About Gene Symbol	




Western blot analysis of TREM2 in various cell lines lysates using TREM2 antibody.




Immunohistochemistry analysis of paraffin-embedded Human tonsil using TREM2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma using TREM2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human gastric cancer using TREM2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse spleen using TREM2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of TREM2 in U251MG using TREM2 antibody, and DAPI (blue).