

DLL3 Rabbit pAb


ABZ1288

Product Overview


Name	DLL3 Rabbit pAb
Accession(Primary)	Q9NYJ7
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:5000-1:20000; IHC-1:50-1:100

Target information (Q9NYJ7)

Gene Name	DLL3
About Gene Symbol	



Western blot analysis of DLL3 in mouse brain and Human fetal brain lysates using DLL3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human cervical cancer using DLL3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.