

CENPM Rabbit pAb

ABZ1289

Product Overview

Name	CENPM Rabbit pAb
Accession(Primary)	Q9NSP4
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:200; IF-1:50-1:200

Target information (Q9NSP4)

Gene Name	CENPM
About Gene Symbol	