


DNA Polymerase mu Rabbit pAb

ABZ1290

Product Overview

Name	DNA Polymerase mu Rabbit pAb
Accession(Primary)	Q9NP87
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:40000

Target information (Q9NP87)

Gene Name POLM**About Gene Symbol**

Western blot analysis of DNA Polymerase mu in K562 lysates using DNA Pol μ antibody.

