


MKNK2 Rabbit pAb

ABZ1291

Product Overview

Name	MKNK2 Rabbit pAb
Accession(Primary)	Q9HBH9
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q9HBH9)

Gene Name MKNK2**About Gene Symbol**

Western blot analysis of MKNK2 in HeLa lysates using MKNK2 antibody.

