

TrpV4 Rabbit pAb

ABZ1292


Product Overview

Name	TrpV4 Rabbit pAb
Accession(Primary)	Q9HBA0
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 98 kDa; Observed MW: 98 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (Q9HBA0)

Gene Name TRPV4

About Gene Symbol



Western blot analysis of TrpV4 in K562 lysates
using TrpV4 antibody.

