


ZBP1 Rabbit pAb

ABZ1296

Product Overview

Name	ZBP1 Rabbit pAb
Accession(Primary)	Q9H171
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:5000-1:20000

Target information (Q9H171)

Gene Name ZBP1**About Gene Symbol**

Western blot analysis of ZBP1 in Raji lysates using
ZBP1 antibody

