

ZBP1 Rabbit pAb

ABZ1297

Product Overview

Name	ZBP1 Rabbit pAb
Accession(Primary)	Q9H171
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:5000-1:20000

Target information (Q9H171)

Gene Name	ZBP1
About Gene Symbol	