

## FTO Rabbit pAb


ABZ1300

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | FTO Rabbit pAb                             |
| <b>Accession(Primary)</b>   | Q9C0B1                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 58 kDa; Observed MW: 58 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; ELISA-1:20000             |

### Target information (Q9C0B1)

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | FTO                                                                                                                                                                                                                                                             |
| <b>About Gene Symbol</b> | FTO (FTO Alpha-Ketoglutarate Dependent Dioxygenase) is a metabolic enzyme involved in translational regulation. It is an important research target in metabolic disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets |



Western blot analysis of FTO in K562 lysates using FTO antibody.