

FOXP3 Rabbit pAb

ABZ1301

Product Overview


Name	FOXP3 Rabbit pAb
Accession(Primary)	Q9BZS1
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Pig
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (Q9BZS1)


Gene Name	FOXP3
About Gene Symbol	




Western blot analysis of FOXP3 in PC-12 lysates using FOXP3 antibody.




Immunohistochemistry analysis of paraffin-embedded Human tonsilla using FOXP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Immunohistochemistry analysis of paraffin-embedded Human bladder cancer using FOXP3 antibody. High-pressure and temperature Tris-EDTA was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human colon using FOXP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human colon using FOXP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human colon using FOXP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.