


TRIM8 Rabbit pAb

ABZ1302

Product Overview

Name	TRIM8 Rabbit pAb
Accession(Primary)	Q9BZR9
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:2000

Target information (Q9BZR9)

Gene Name TRIM8**About Gene Symbol**

Western blot analysis of TRIM8 in Human placenta and Human fetal brain lysates using TRIM8 antibody.

