

TPSD1 Rabbit pAb

ABZ1303

Product Overview

Name	TPSD1 Rabbit pAb
Accession(Primary)	Q9BZJ3
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 27 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (Q9BZJ3)

Gene Name TPSD1

About Gene Symbol



